LUBRAGEL® - Sterile Lubricant Gel with Lidocaine-

INSTRUCTIONS FOR USE - Please read carefully.

DESCRIPTION

LUBRAGEL is a sterile, single-use syringe or accordion tube, pre-filled with a water-soluble gel containing Lidocaine and Chlorhexidine.

INGREDIENTS

100g of LUBRAGEL contains:
0.05g Chlorhexidine Digluconate, 2g Lidocaine Hydrochloride, Purified Water, Propylene Glycol, Hydroxyethyl cellulose, Methyl Hydroxybenzoate, Propyl Hydroxybenzoate.

INTENDED PURPOSE

LUBRAGEL is a medical device that helps to prevent trauma being caused to the patient during catheterisation procedures or other urethral procedures (e.g. endoscopy and cystoscopy) by effective lubrication. LUBRAGEL contains Chlorhexidine as an antiseptic to help reduce the risk of infection and Lidocaine as a local anaesthetic to help reduce pain. By the use of LUBRAGEL the patient is relaxed and iatrogenic injuries due to spasticity or discomfort are minimized.

INDICATION FOR USE

LUBRAGEL is designed to lubricate the urethra prior to insertion of a urethral catheter and other urological medical devices including cystoscopes. It can also be used as a lubricant gel for rectal and colonic applications. It can also be used as a lubricant for all forms of endoscopy and to prevent injuries to the colon and rectum. The transparent gel leaves the optics of the instruments visibly clear.

CONTRAINDICATIONS

The gel must not be used in patients with known hypersensitivity to any of the ingredients. LUBRAGEL must not be applied when other medical devices or medicines containing local anaesthetic agents are used at the same time. It should not be used in patients who have damaged or bleeding mucous membranes because of the risk of systemic absorption of Chlorhexidine and Lidocaine. Do not use in children below 2 years.

WARNINGS AND PRECAUTIONS

LUBRAGEL must be used only under the supervision or according the instruction of health care professionals.

LUBRAGEL is not suitable if a patient

- ever had a reaction to a local anesthetic,
- is allergic or hypersensitive to Parabens, Chlorhexidine Digluconate or any of the other ingredients,
- has a damaged or bleeding urethra.

Care should be taken when using LUBRAGEL with patients

- who have heart problems or are taking medication for treating irregular heartbeat,
- who have liver problems or
- who are epileptic.

The patient might feel a slight stinging just after the gel is used. This stops as soon as the anesthetic starts to work. If patients experience any reaction to the lubricant gel, they should inform their doctor as soon as possible. If any of the side effects gets serious or if any side effects not listed in this leaflet occur, the use should be stopped and the doctor, pharmacist or manufacturer should be informed.

The gel must not be used orally. If the gel has been used orally, care should be taken when chewing or swallowing. Numbness of the tongue or mouth are can lead to a bite injury. Contact with the eyes must be avoided because numbness in the eyes may prevent the patient from noticing if they got foreign bodies in their eye.

LUBRAGEL is not to be used if the packaging is damaged. Do not use after expiry date. LUBRAGEL is for single use only. Not to be used for i/v and i/m injections. Keep out of reach of children.

PREGNANCY AND LACTATION

Before using LUBRAGEL, patients should always be asked if they are pregnant. It should only be administered during pregnancy or breastfeeding on the advice of a doctor.

EFFECTS ON ABILITY TO DRIVE AND USE MACHINES

The ability to drive and operate machinery may be slight impaired after the use of LUBRAGEL. If affected, patients should be advised not to drive or use machinery.

SIDE EFFECTS

Undesirable effects of Lidocaine are possible where there is severe injury to the mucosa and absorption may occur. Local hypersensitivity reactions may occur in rare cases such as redness, stinging, itching or blistering and/or systemic reactions to Lidocaine and/or Chlorhexidine. There is also a risk of severe reactions including a drop in blood pressure, dizziness, nausea, shortness of breath, bradycardia, convulsions and anaphylactic shock.

DRUG INTERACTIONS

Depending on the absorption of Lidocaine, interactions can be seen when used with the following medications:

- Propranolol: Reduction in plasma clearance of Lidocaine,
- Cimetidine: Reduction in plasma clearance of Lidocaine,
- Antiarrhythmic products: Increase in the toxicity of Lidocaine,
- Phenytion or barbiturates: Reduction in plasma levels of Lidocaine.

Specified interactions can be seen in the long-term use and repeated high doses.

When administered as recommended, there are no reported clinically significant interactions.

DOSAGE RECOMMENDATIONS

Adults; Max. 800 mg Lidocaine during a 24 hour period, i.e. three 11-ml-syringes or six 6-ml- syringes or four 12.5 g accordion tubes.

Children (between 2 and 15 years): Max. 0,3 ml gel/kg B.W (ca. 6mg Lidocaine/kg) is recommended per procedure and no more than four doses should be administered within 24 hours.

Cystoscopy: The entire urethra including the sphincter externus must be covered with LUBRAGEL and locally anaesthetised; another syringe may be required for this. The decision about the amount of LUBRAGEL to be used is taken by the physician.

OVERDOSE

LUBRAGEL must not be applied when other medical devices or medicines containing local anaesthetic agents are used at the same time.

Symptoms: In the event of excessive absorption of Lidocaine into the bloodstream, symptoms may include central nervous system effects and cardiovascular reactions. If signs of an overdose occur, a doctor is to be consulted immediately. For further questions on excessive use or use outside of these instructions, a medical doctor or pharmacist should be contacted.

Treatment: The treatment of a patient suffering from systemic toxicity of Lidocaine consists of arresting the convulsions and ensuring adequate ventilation with oxygen, if necessary by assisted or controlled ventilation (respiration).

HOW TO USE

It is recommended that this product is used as part of an aseptic technique.

- Cleanse the opening to the urethra and surrounding area prior to use.
- Evenly peel back the paper backing and remove the sterile syringe or accordion tube.
- Syringe: Before removing the cap from the end of the syringe, loosen the plunger by gently pressing it. Then remove the cap.
- Accordion tube: Bend the tip, if possible still in the blister pack. Remove the tip completely to avoid accidental insertion into the urethra.
- Apply a drop of gel to the urethral opening to make initial insertion easier;
- Insert the nozzle into the urethral opening and press the plunger of the syringe slowly to release the necessary amount of the gel.

NOTE: LUBRAGEL should be applied to the desired body area, not directly onto the device.

The lubricating characteristics start to take effect at the same time of application. The onset of the anesthetic effect takes 3-5 minutes. The duration of anaesthesia is approximately 20-30 minutes.

STORAGE / SHELF LIFE / DISPOSAL

Store at 5-30°C until expiry date. Keep away from sunlight and keep dry. The shelf life of the product is 3 years. The expiration date is printed on the blister packaging and on the box. Any remaining gel in the syringe/accordion tube must be disposed according to local policies and procedures.

STERILIZATION

LUBRAGEL is sterilized using irradiation after the packaging process.

PACKAGING

LUBRAGEL is available in syringes prefilled with 6ml (ca. 6g) or 11ml (ca. 11g) gel and accordion tubes filled with 12,5g gel. The scale on the syringe does not have a measurement function; it is only for the orientation of the user.

LG-L 006 LUBRAGEL	6 ml (25x6 ml/Box)
LG-L 011 LUBRAGEL	11 ml (25x11ml/Box)
LG-L 12.5 LUBRAGEL	12,5g (25x12,5g/Box)

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4M Medical GmbH