UROLUBE® - Sterile Lubricant Gel with Lidocaine - Paraben-free and Chlorhexidine-free

INSTRUCTION FOR USE - Please read carefully.

DESCRIPTION

UROLUBE is a sterile, single-use syringe or accordion tube, pre-filled with a water-soluble gel containing Lidocaine. It is Paraben-free and free of Chlorhexidine-free.

INGREDIENTS

100g of UROLUBE contains:

Purified Water, Propylene Glycol, Hydroxyethylcellulose (Lubricative) and 2g Lidocaine Hydrochloride (Local anaesthetic).

INTENDED PURPOSE

UROLUBE is a medical device that helps to prevent trauma being caused to the patient during catheterization procedures or other urethral procedures (e.g. endoscopy and cystoscopy) by effective lubrication. UROLUBE contains Lidocaine as a local anaesthetic to help reduce pain. By the use of UROLUBE the patient is relaxed and iatrogenic injuries due to spasticity or discomfort are minimized.

INDICATION FOR USE

UROLUBE is designed to lubricate the urethra prior to insertion of a urethral catheter and other urological medical devices including cystoscopes. It can also be used as a lubricant gel for rectal and colonic applications. It can also be used as a lubricant for all forms of endoscopy and to prevent injuries to the colon and rectum. The transparent gel leaves the optics of the instruments visibly clear.

CONTRAINDICATIONS

UROLUBE must not be used in patients with known hypersensitivity to the active ingredients or any of the excipients. It should not be used in patients who have damaged or bleeding mucous membranes because of the risk of systemic absorption of the lidocaine hydrochloride. **Do not use in children below 2 years.**

WARNINGS

UROLUBE must be used only under the supervision of or according to the instruction of health care professionals.

UROLUBE is not suitable if a patient

- ever had a reaction to a local anaesthetic,
- is allergic or hypersensitive to any of the other ingredients,
- has a damaged or bleeding urethra or membrane.

Care should be taken when using UROLUBE with patients

- who have heart problems or are taking medication for treating irregular heartbeat,
- who have liver problems,
- who are epileptic,
- who are pregnant or breast feeding.

The patient might feel a slight stinging just after the gel is used. This stops as soon as the anesthetic starts to work. If patients experience any reaction to the lubricant gel, they should inform their doctor as soon as possible. If any of the side effects gets serious or if any side effects not listed in this leaflet occur, the use should be stopped and the doctor, pharmacist or manufacturer should be informed.

UROLUBE must not be used orally. If the gel has been used orally, care should be taken when chewing or swallowing. Numbness of the tongue or mouth can lead to a bite injury. Contact with the eyes must be avoided because numbness in the eyes may prevent the patient from noticing if they got foreign bodies in their eye.

UROLUBE is not to be used if the packaging is damaged. Do not use after expiry date. UROLUBE is for single use only. Not to be used for i/v and i/m injections. Keep out of reach of children.

PREGNANCY AND LACTATION

Before using UROLUBE, patients should always be asked if they are pregnant. It should only be administered during pregnancy or breastfeeding on the advice of a doctor.

EFFECTS ON ABILITY TO DRIVE AND USE MACHINES

The ability to drive and operate machinery may be slight impaired after the use of UROLUBE. If affected, patients should be advised not to drive or use machinery.

SIDE EFFECTS

In spite of the proven wide safety range of UROLUBE undesirable effects of the local anaesthetic, lidocaine is possible where there is severe injury to the mucosa and absorption may occur. Examples are anaphylaxis, fall in blood pressure, bradycardia or convulsions.

DRUG INTERACTIONS

Depending on the absorption of Lidocaine, these interactions can be seen when used with the following medications:

Propranolol: Reduction in plasma clearance of lidocaine,
Cimetidine: Reduction in plasma clearance of lidocaine,
Antiarrhythmic products: Increase in the toxicity of lidocaine,
Phenytoin or barbiturates: Reduction in plasma levels of lidocaine.
Specified interactions can be seen in the long-term use and repeated high doses.
When administered as recommended, there are no clinically significant interactions reported.

DOSAGE RECOMMENDATIONS

Adults; Max. 800 mg Lidocaine during a 24-hour period, i.e. three 11-ml-syringes or six 6-ml-syringes or four 12,5 g accordion tubes.

Children (between 2 and 15 years): Max. 0,3 ml gel/kg B.W (ca. 6mg Lidocaine/kg) is recommended per procedure and no more than four doses should be administered within 24 hours.

Cystoscopy: The entire urethra including the sphincter externus must be covered with UROLUBE and locally anaesthetised; another syringe may be required for this.
The decision about the amount of UROLUBE to be used is taken by the physician.

OVERDOSE

UROLUBE must not be applied when other medical devices or medicines containing local anaesthetic agents are used at the same time.

Symptoms: In the event of excessive absorption of Lidocaine into the bloodstream, symptoms may include central nervous system effects and cardiovascular reactions. If signs of an overdose occur, a doctor is to be consulted immediately. For further questions on excessive use or use outside of these instructions, a medical doctor or pharmacist should be contacted.

Treatment: The treatment of a patient suffering from systemic toxicity of Lidocaine consists of arresting the convulsions and ensuring adequate ventilation with oxygen, if necessary by assisted or controlled ventilation (respiration).

HOW TO USE

It is recommended that this product is used as part of an aseptic technique.

- Cleanse the opening to the urethra and surrounding area prior to use.
- Evenly peel back the paper backing and remove the sterile syringe or accordion tube.
- Syringe: Before removing the cap from the end of the syringe, loosen the plunger by gently pressing it. Then remove the cap.
- Accordion tube: Bend the tip, if possible still in the blister pack. Remove the tip completely to avoid accidental insertion into the urethra.
- Apply a drop of gel to the urethral opening to make initial insertion easier;
- Insert the nozzle into the urethral opening and press the plunger of the syringe slowly to release the necessary amount of the gel.

NOTE: UROLUBE should be applied to the desired body area, not directly onto the device.

The lubricating characteristics start to take effect at the same time of application. The onset of the anaesthetic effect takes 3-5 minutes. The duration of anaesthesia is approximately 20-30 minutes.

STORAGE / SHELF LIFE

Store at 5-30°C until expiry date. Keep away from sunlight and keep dry. The shelf life of the product is 3 years. Do not use after expiry date. The expiration date is printed on the blister packaging and on the box. Any remaining gel in the syringe/accordion tube must be disposed according to local policies and procedures.

STERILIZATION

UROLUBE is sterilized using irradiation after the packaging process.

PACKAGING

UROLUBE is available in syringes pre-filled with 6ml (ca. 6g) or 11ml (ca. 11g) gel and accordion tubes filled with 12,5g gel. The scale on the syringe does not have a measurement function; it is only for the orientation of the user.

UL-L 006 UROLUBE	6 ml (25x6 ml/Box)
UL-L 011 UROLUBE	11 ml (25x11ml/Box)
UL-L 12.5 UROLUBE	12,5 g (25x12,5g/Box)

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4M Medical GmbH